

TANYX

A NEW CONCEPT IN PAIN TREATMENT

Effective against:

Muscle Aches

Osteoarthritis

Menstrual Cramps

Back Pain

Tendinitis



Portable & Easy to use



Drug Free



No Side Effects



Organic Action



Provides 30 Sessions

QUICK MANUAL

01

Clean dirt and oils from application area of the skin and allow to dry. If necessary, trim excessive hair for better adhesion.

02

Before turning the device ON remove and discard the **Electrode Gel Pad Sticker** (). Stick one gel pad onto each electrode of TANYX® device. *This Electrode Gel Pad will remain until battery is depleted.*

03

Remove and save the **Reusable Gel Pad Sticker (Protection)***. Apply TANYX® to the desired area ensuring device adheres to skin by pressing firmly on both sides. **This protective sticker will be used to store the device when not in use.*

04

3 SEC.

ON: Press and hold the button until the indicator light starts blinking and choose the intensity: L (Low), M (Medium), H (High) by pressing the respective button.
OFF: Press and hold the button until the indicator light stops blinking.

05


Before removing TANYX® turn device, hold down the skin with one hand while gently pulling the device up with the other hand. *Ensure Electrode Gel Pad does not detach from the device.*

06

Apply **Reusable Gel Pad Sticker (Protection)** you set aside onto each gel pad. Place TANYX® device inside the **blister pack** provided and store it for future uses. *Never detach Electrode Gel Pads from TANYX® device.*

RECOMMENDED APPLICATIONS AREAS

Front



Back

